

Medela Clario Toni suction pump

Pump Specifications	
Item numbers	AC/DC version
	014.0022 100–240V/50–60Hz AUS/NZ plug A
	014.0023 100–240V/50–60Hz Euro plug A
	014.0024 100–240V/50–60Hz Euro plug B
	014.0025 100–240V/50–60Hz UK Plug A
	Deluxe version (AC/DC + SoftCase)
	014.0026 100–240V/50–60Hz AUS/NZ plug A
	014.0027 100–240V/50–60Hz Euro plug A
	014.0028 100–240V/50–60Hz Euro plug B
	014.0029 100–240V/50–60Hz UK Plug A
Instructions for use	200.8400 EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A) 200.8401 ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)
Service manual	200.4001 EN, DE, FR
Max. vacuum [kPa/mmHg] Measured at 0 m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).	Tolerance: min. –9/–70 ± 2.6/20 med. –15/–115 ± 4/30 max. –80/–600 – 12/90
Flow [l/min]	15 l/min. ± 10%
Dimensions [mm]	223 x 255 x 95
Weight [kg]	2.0
Drive unit	Membrane
Power supply	100–240V, 50–60 Hz
Power consumption [V]	12V, 3.4 A
Noise level [dB(A)]; measured in Medela acoustic room at 1 m distance	63.2
Protection against electric shock	Type BF
Medical Device Classification – According to 93/42/EEC	Class IIa
Protection against electric shock	Class II Equipment
Ingress Protection	IP21
Battery	NiMH
Operation from battery (at full charge)	>50 minutes
Charging time of battery	5 hours
Overflow protection	Threefold hygienic design: – Floater in canister prevents overflow – SafetyChamber – Membrane system No overflow can occur into the pump
Warranty [years]	2; Battery 6 months
Operating temperature, relative humidity, pressure	5 to 40 °C, 15–93 %, 70–106 kPa
Transport/storage temperature, relative humidity, pressure	–20 to +50 °C, 15–93 %, 70–106 kPa
Made in	Switzerland


Accessories	Item number
Canister set	074.0001
SafetyChamber Clario Toni	074.0030
Disposable Fingertip tubing	074.0003
Car connection cable	074.0017
SoftCase carrying bag	074.0007
Vacuum gauge	074.0006

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.



Medical Vacuum Technology
for Healthcare Professionals

Please contact us or your local Medela
representative for details.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medelahealthcare.com

Local contact: